

** Note: Orders are filled within 24-48 hours from time of receipt and delivered Monday - Friday**

For Supply Questions Call: (503) 814-5227, opt. 1 Fax Complete Orders to: (503) 480-4853

SUPPLY REQUEST FORM

dered by:	Order Date: Date Filled:
Salem Health Regional Laboratory Supplies	Aurora/Pacific Pathology Supplies
OOD COLLECTION SUPPLIES	CYTOLOGY SUPPLIES
Vacutainer, 2.7 ml blue top (100 tray)	Pap kit, Sure Path [™] (1tray/25 ea.)
Vacutainer, 4.5 ml green top (100 tray)	Sure Path™ Brushes/Spatula (1 pk/25 ea
	Brooms, Blue Sure Path™ (1 pk/25 ea.)
Vacutainer, 4.0 ml lavender top (100 tray) Vacutainer, 5 ml gold top (100 tray)	
(Serum separator tube, no anticoagulant)	Pap kit, Thin Prep® (1 tray/25 ea.)
Vacutainer, 6 ml red top (100 tray)	Thin Prep® Brushes/Spatulas (1 pk/25
(No anticoagulant)	Brooms, Purple Thin Prep® (1pk/25 ea.
Vacutainer needles, 21 gauge (100 box)	
Vacutainer Needles, 22 gauge (100 box)	Slides, (frosted on one end) (box)
Needle holders, disposable (25 per bag)	Slide mailer, 5-slide plastic (ea.)
	Tubes w/ CytoLyte solution 50ml (ea.)
CROBIOLOGY SUPPLIES	Pap slide holders w/alcohol (ea.)
	RPMI (must be refrigerated) (ea.)
Blood culture bottle (Aerobic)	Afirma Thyroid FNA (ea.)
Blood culture bottle (Anaerobic)	
Blood culture bottle (Pediatric)	PATHOLOGY SUPPLIES
Blood Culture collection Insert	1 oz. formalin in 30 ml container (1 tray/36 ea.)
(for use when collecting tubes after blood culture)	2 oz. formalin in 60 ml container (1 tray/25 ea.)
Aerobic transport (E Swab) (EACH)	4 oz. formalin in 120 ml container (1 tray/25 ea.
BD Affirm (Box of 10)	Formalin bottles, 1 gallon
	Formalin, 15 ml Container (Box)
Universal Transport media (herpes)	Specimen container, 8 ounce
	Specimen container, 16 ounce
RSV & Influenza A/B Collection Kits	Specimen container, 32 ounce
Pertussis Collection Kits	Specimen container, 1/2 gallon
	Specimen container, 1 gallon
Stool collection kits with instructions	4 oz. Sterile Cup (ea.)
	Aurora/PPA Specimen Transport Bags
	Quart (pk.) Gallon (pk.)
NERAL SUPPLIES	FORMS
	Aurora/Pacific Pathology:
Urine Kit - Sterile cup, C&S w/ straw, UA pres. tube	Gynecological Requisitions
(kits come in bags of 10 kits ea.)	Surgical Requisitions
4 oz. Urine container, sterile (no preservative) (bags of 25 ea.)	
BD Urinalysis Preservation Tubes (ea.) (yellow and red top)	Salem Health:
BD Urine Transfer Straw w/ C&S tube (ea.) (gray tube)	Salem Health Laboratory Test Requisitions
24 Hr. Urine Container (ea.)	Supply Request Form
Gonorrhea/Chlamydia collection kit (CT/NG Swab)	NextGen Printer Paper
Sonormea/Chiamydia Conection Kit (CT/NG 3wab)	OTHER
Salem Hospital Specimen Transport Bags (pkg of 100)	VIIIEN
Regular Specimen (Quart) Regular Specimen (Quart) Large (Gallo	on)
Large (Gallo	<u> </u>
end completed form with courier or fax to : 503-480-1853	
rovision of Supplies and Equipment	

It is our policy to provide supplies and equipment directly related to the services we provide for our clients. Services and supplies to be used by clients for other purposes will either be billed at fair market value or we will provide a referral to a high quality independent vendor.